



3DSlicer user guidelines

Step 1 - Download 3DSlicer **4.13.0 version** from google.


Step 2 - Open the application by double clicking on the **3DSlicer icon**.




Step 3 - Click on the **View icon (Green box)** and choose the **Extension manager** in the dropdown menu.




Step 4 - Click on the **Install Extensions(Greenbox)**




Step 5- In the search bar type **Qualitative reporting** and Install it.



Step 6 - Drag and drop the Dicom and its corresponding Seg files (Green arrow) to the **3D Slicer** (Red arrow).



Step-7 - Select both the uploaded files and click **load** (**Bottom green box**).



Step 8 - Annotations can be viewed in the **3DSlicer screen**.

